

## Checklist for personalized NIBS for L2 learning - Electrodes placement

- **Montage**

- ☐ FC3: Anode (Red) & T7: Cathode (Blue)
- ☐ Stim. elec. Vicinity: FFC3h, FFC5h, FCC3h, FCC5h & None

- **Preparation tACS Electrodes**

- ☐ **EEG cap:** Measure brain size (circumference)  $\Rightarrow$  Choose and put on EEG cap  $\Rightarrow$  Adjust position (10-20 system)  $\Rightarrow$  **Fasten the chin strap**
- ☐ **Location marking:** Use EOG sticker to mark FC3 & T7  $\Rightarrow$  Remove the caps  $\Rightarrow$  Remove the hair in marked area
- ☐ **Stimulator:** OFF  $\Rightarrow$  Connect the wire  $\Rightarrow$  Apply Ten20 gel on the tACS electrodes  $\Rightarrow$  Place it in the marked area  $\Rightarrow$  Turn ON stimulator
- ☐ **Impedance check (10k Ohm):** Add more paste or press until the impedance is under  $10k\Omega$  and carefully pressing on the tACS electrodes to avoid tACS electrodes sliding.
- ☐ **Remove the excess gel with a cotton swab**
- ☐ Turn OFF the stimulator

- **Mounting the EEG cap**

- ☐ Carefully and gently put on the cap and MAKE SURE the tACS electrodes are not shifted and no Ten20 gel leaked out
- ☐ Fasten the strap of the EEG cap

- **Preparation for EEG electrodes**

- ☐ **Gel electrode:** ground&ref  $\Rightarrow$  Vicinity of the tACS electrodes  $\Rightarrow$  Rest
- ☐ **Note for vicinity:** Inject the gel with the needle tip pointing in a direction away from the tACS electrode and apply as LESS GEL as possible.
- ☐ First round
- ☐ Second round: In the second round use the **cotton swab for the vicinity**
- ☐ **EEG impedance check:** gel until green/yellow
- ☐ **!!! Hardware check:** DC mode + 1Kz sampling frequency
- ☐ **!!! Stimulator setting:** Signal out ON, REMOTE ON
- ☐ **Saturation test,** i.e., apply 1mA sinusoidal signal and observe any saturation of EEG electrodes from recording signal. (safe\_start.py & safe\_end.py)
- ☐ **Check EEG AND tACS impedance again and start recording.**

